



	<b>Experiment title:</b> High-pressure XANES study of sulphur speciation in crude oil.	<b>Experiment number:</b> ES-181
<b>Beamline:</b> ID12	<b>Date of experiment:</b> from: 12.02.2014 to: 17.02.2015	<b>Date of report:</b> 25.02.2015  <i>Received at ESRF:</i>
<b>Shifts:</b> 18	<b>Local contact(s):</b> Fabrice Wilhelm	

**Names and affiliations of applicants (\* indicates experimentalists):**  
**Innokenty Kantor (ESRF)**  
**Leonid Dubrovinsky (University of Bayreuth, Germany)**


## Report:

Experiments under high pressures at a low energy range (~2.5 keV) imposes a lot of technical challenges, and most of them were successfully resolved. Ultra-thin (30 microns) diamond windows demonstrate sufficient strength to maintain pressure at least up to 7 GPa (in one test run), while still having enough X-ray transmittance.


However, two main technical problems are still not resolved completely: sulphur contamination of the sample environment and radiation damage of the soft matter samples.

As was found experimentally, the commercially available diamond windows contain a detectable amount of sulphur in a sulfate form (Fig. 1), mostly at the laser-cut edges. This contamination occurs due to the windows cleaning in a sulphuric acid (that penetrates inside the microcracks) during manufacturing, and it cannot be easily removed. To reduce the contamination effect due to the X-ray beam tails we used CRL focusing instead of a mirror, which helped to a certain extent, but at the same time also enhanced the radiation damage of the sample due to the harmonics (see Fig. 1).

Nevertheless, we demonstrated a possibility to obtain a high-quality XAS spectra in a tender X-ray domain under high pressure, also from natural samples with a low (~1.5%) concentration of sulphur (Fig. 2). The observed changes of S K-edge XANES spectra in crude oil as a function of pressure were very subtle, and these small changes could be probably attributed to the radiation damage of the sample.



Left (Fig. 1). The reference S K-edge XANES spectrum of the crude oil from West Siberia (top blue curve) and a sequence of spectra of the same sample inside the diamond anvil cell taken as a function of time (top to bottom). The star mark shows a feature at  $\sim 2482$  eV, coming from the sulphur contamination of the diamond windows. Two arrows indicate spectral features, specific to the sample radiation damage: a growing pre-edge peak at  $\sim 2468$  eV and a reducing feature at  $\sim 2454$  eV.



Right (Fig. 2). XANES spectra of the natural oil sample as a function of pressure at ambient temperatures. No changes in the sulphur oxidation state were observed, and the small spectra changes are probably attributed to the radiation damage.