

Anyhow, we proved the feasibility of the experiment in the proposed experimental condition as we managed to collect XANES spectra of $\text{Fe}(\text{TPA})(\text{OTf})_2$ in 10 ms before mixing it with H_2O_2 . As the reaction we want to study occurs in 100 ms we are confident that we will manage to characterize the nature of the intermediates that are formed. We managed to improve the quality of the spectra by collecting a spectrum of the cell windows that was then subtracted from the sample spectrum. The results are shown in Figure 1 where we show the raw Fe K-edge XANES spectrum of a 20 mM solution of $\text{Fe}(\text{TPA})(\text{OTf})_2$ in acetonitrile together with the absorption spectrum of the empty cell.




Figure 1. In the upper panel we show the raw XANES spectrum at the Fe K-edge of a 20 mM solution of $\text{Fe}(\text{TPA})(\text{OTf})_2$ in acetonitrile, collected in 10 ms (black line) together with the spectrum of the empty cell collected with the same statistics. In the lower panel we show the XANES spectrum obtained by subtracting the two spectra showed above.

The noise of both spectra is obviously due to the strong absorption of the windows that is also responsible for the distortion in the interval between 7150 and 7250 eV. In the lower panel of the figure we show the XANES spectrum obtained after subtraction, and the signal to noise ratio is very much improved. Moreover, the edge position is perfectly detected thus allowing us to determine variations of the oxidation state of Fe that is one of the key points

of this project. Also the XANES oscillations are not distorted anymore and they can be analysed on a quantitative ground.

We believe that the development of this experimental technique has the potential of a much wider application as a novel mechanistic tool to investigate metal-catalyzed redox reactions. Indeed, there is currently a lively debate on the nature of the actual oxidizing species in several oxidation catalysts (mostly on iron based ones), and from a broader perspective high oxidation states of transition metals are receiving increasingly interest for their unexpected and promising reactivity. This technique would also provide a complementary tool to characterize the intermediates involved in catalytic systems where conventional spectroscopic methods (Mossbauer, EPR) cannot work (such as in manganese catalyzed oxidations, as Mn is EPR and Mossbauer silent).