



Experiment title: Structural Studies of the Restriction
Endonuclease Bgl1.

**Experiment
number:
LS668**

Beamline:
D2AM,
BM14

Date of experiment:
from: 9th June 1997 to: 10th June 1997
from: 11th June 1997 to: 13th June 1997

Date of report:
13th August

Shifts:
9

Local contact(s):
Eric Fanchon, Andy Thompson

Received at ESRF:
22 AOUT 1997

Names and affiliations of applicants (* indicates experimentalists):

Phillips S.E.V., Newman M.*, Smith L.*, Vincombe E.*, Wilmot C.M.*
Department of Biochemistry and Molecular Biology,
University of Leeds,
Leeds
LS2 9JT
United Kingdom.

Report:

We collected four complete Bgl1/DNA complex data sets during our beam time on D2AM three of which are useful and are being used in ongoing structural determinations:

Native and PCMBS derivative data sets were collected for a Bgl1/17mer complex at wavelengths of 1.00911 and 1.00465Å respectively. The derivative data has not been processed but the native data are good to 2.24Å In addition, 2.24Å native and PCMBS derivative data sets were collected for a Bgl1/18mer complex at wavelengths of 1.00917 and 1.00465Å respectively.

In total, we collected seven complete Bgl1/DNA complex data sets during our beam time on BM14 six of which are useful and are being used in ongoing structural determinations.

Native, PCMBS and hexachloro osmate derivative data sets were collected for the Bg11/17mer complex at wavelengths of 0.919, 0.919 and 1.14Å. The 2.2Å native data have an R_{symm} of 0.055 for all data (0.078 in the last resolution shell), the 2.3Å PCMBS data have an R_{symm} of 0.039 for all data (0.107 in the last resolution shell) and the 2.8Å hexachloro osmate data have an R_{symm} of 0.031 for all data (0.054 in the last resolution shell).

2.3Å data sets were also collected for three brominated Bg11/17mer complexes at a wavelength of 0.919Å. These data have R_{symm} 0.034, 0.038, 0.033 for all data (0.065, 0.081, 0.067 in the last resolution shell respectively).

A hexachloro osmate derivative data set was collected for the Bg11/18mer complex at a wavelength of 1.14Å, but these data have not been processed.

The useful data sets have been incorporated into the MIRAS phase calculation using the program SHARP, resulting in a superb electron density map for the Bg11/17mer complex. A model has been built, including the DNA and all but one residue of the protein and is currently being refined.

We are pleased with the data collected during this beam time and would like to thank the beam line staff, in particular Eric Fanchon and Andy Thompson, for their assistance.